



# Chemistry Behind

## BEAUTIFIL Injectable

**BEAUTIFIL Injectable X & XSL** feature a novel filler technology that involves controlled miniaturization and uniform nano-silanization for higher filler load and optimal filler-resin ratio for superior strength and durability.


The proprietary bioactive nano S-PRG (Surface Pre-Reacted Glass ionomer) filler technology provides an optimal balance of light transmission / diffusion for a perfect shade match and bioactive benefits. The homogenous nano-filler complex facilitates an easy and lasting self-polishing natural gloss.




Photographs courtesy of Prof. Koji Yamamoto of the ASAHI University

## Exceptional Aesthetics that Last


**BEAUTIFIL Injectable X** exhibits a high level of natural aesthetics with life-like shades that merge seamlessly with natural tooth for invisible restorations. Shades are stable with no shift in translucency pre-and post-cure.




### Light transmission and diffusion




### BEAUTIFIL Injectable X (15 shades)




B1 B2 C2 D2



### BEAUTIFIL Injectable XSL (13 shades)




B1 B2 A20 A30



## 3 Levels of Opacity to meet all clinical needs

- Opaque shades to mimic deep dentin opacity
- Universal shades for universal application
- Universal enamel to simulate translucency of natural enamel

The shades A1-D2 correspond to the VITA® Classical Shades.



**BEAUTIFIL Injectable X & XSL** exhibit effortless polish to an unparalleled natural luster with self-polishing ability for lasting aesthetics with higher stain resistance.

High gloss observed post-cure  
(without polishing)