

## The 2<sup>nd</sup> Generation Giomer Aesthetic Dental Restorative Material

BEAUTIFIL II



**BEAUTIFIL II** the 2<sup>nd</sup> generation Giomer light-cure, fluoride releasing direct aesthetic restorative material was developed using S-PRG Technology (Surface Pre-Reacted Glassionomer) to provide superior aesthetics with excellent shade match, enhanced handling characteristics and anti-plaque effect while maintaining strength and durability.

Confidence starts here...



Photographs courtesy of Dr. Sushil Koirala

### Created for ...

Direct restorations that require optimum aesthetics and biocompatibility such as

- Direct cosmetic tooth modifications
- Restoration of Class III, IV, and V cavities
- Restoration of Cervical erosion and root caries
- Restoration of Class I and selective Class II cavities
- Repair of fractured incisal edges
- Laminate veneers and core build-up

### Predictable Aesthetics



Create Indistinguishable restorations with optical characteristics of natural teeth...

The **BEAUTIFIL II** filler structure has been developed to simulate the internal structure of natural teeth with ideal light transmission and optical characteristics.

The moderate translucency and light transmission of enamel combined with the light-diffusion of dentin offers predictable aesthetics with a close shade match to natural teeth.



Photograph courtesy of Dr. Sushil Koirala



### Versatility...

Excellent natural shade reproduction can be achieved with a chameleon effect, using a single shade that blends well with surrounding teeth making the restoration undetectable. In aesthetically demanding cases additional shades can be used to achieve exceptional results.

**BEAUTIFIL II** is available in 14 shades –A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, Inc, BW and Opacious Dentin shades (A2O, A3O)

### Recommended Use of Shades

**Universal shades – A1, A2, A3, A3.5, A4, B1, B2, C2**

Designed specially for the single shade restorative technique

**Opacious Dentin shades – A2O, A3O**

Specially created for restoration with layering technique using multiple shades


**Enamel shade-Inc**

Unique Incisal shade with high translucency for aesthetically demanding cases

**Bleach White shade – BW**

Ideal for restoring bleached teeth

### Fluorescence close to natural teeth




### Exceptional Radiopacity



1.7 times of Enamel and 3 times higher than Dentin

### Shade Stability



Negligible difference in shade before and after curing


### Optimum Polishability

In addition to the S-PRG fillers discrete nano fillers (10~20nm) have been included in the filler structure of **BEAUTIFIL II** to obtain a filler load of 83.3 wt% for fast and easy polishing with an outstanding surface lustre which retains over time.

### Outstanding Handling Characteristics

Excellent handling has been achieved with the inclusion of nano particles making **BEAUTIFIL II** a user-friendly material that is easy to sculpt with no slump, non-sticky and wets preparation well with good thixotropic qualities.

The "One Touch Cap" of the syringe has been enhanced for easy dispensing and delivery during build-up of restoration.



**BEAUTIFIL II**  
REFILLS: Shades:  
A1, A2, A3, A3.5,  
A4, B1, B2, B3, C2,  
C3, Inc, BW, A2O, A3O

each 4.5 gm  
\$49.50



**BEAUTIFIL II**  
6 COLOUR SET 6 x Beautifil II  
Syringe, 4.5gm A2, A3, A3.5,  
B2, A3O, Inc, Dura-White  
Stone, Super-Snap Singles,  
OneGloss Trial Set & Paper Pad

\$247.50 / set



**Buy 5 get**  
**1 FREE**  
(same item)