



ZiLMaster

The Key to Polishing Modern Ceramics


ZiLMaster is a colour co-ordinated, three-step diamond impregnated polisher, specially designed with a uniform distribution of diamond particles in a resilient silicone binder to provide the smoothest surface finish and high gloss polish on modern all-ceramic restorations such as **Zirconia** and **Lithium Disilicate**.

ZiLMaster comes in a variety of shapes in 3 specific grits for both intra-oral and chairside application.


CONTOURING WITH DURA-GREEN DIA




Contouring of oxide-ceramics such as Zirconia involves shaping as well as the removal of larger surface scratches. Dura-Green DIA is a specially developed vitrified diamond impregnated abrasive that yields a smoother surface with virtually no chipping and less heat generation.



Surface image (x 400) after contouring with Dura-Green DIA and Diamond Point Fine





Contoured with Dura-Green DIA



Contoured with Diamond Point Fine

PREDICTABLE STEPS TO A LASTING GLOSS

Material	Extra-oral / Intra-oral	Recommended procedure from contouring to polishing			
		Contouring	Finishing	Polishing	Super Polishing
Zirconia	Extra-oral Intra-oral	Dura-Green DIA	ZiLMaster Coarse	ZiLMaster Medium	ZiLMaster Fine
Lithium Disilicate (or other ceramics)	Extra-oral Intra-oral		*Only in delicate areas such as margins or edges of inlays/onlays/laminate veneers.	ZiLMaster Medium	



(Internal data: Shofu R&D)