


## BIOACTIVE BENEFITS

### All round protection ...

Proprietary bioactive S-PRG Fillers incorporated in **BEAUTIFIL II LS** exhibit acid-neutralizing ability to maintain an acid neutral and fluoride rich environment. Research has shown that Giomer based restorative materials including **BEAUTIFIL II LS** have an anti-plaque effect on the restoration surface for long-term protection and to prevent potential gingival irritation.



Photographs courtesy of Prof. Koji Yamamoto of the ASAHI University

Photographs courtesy of Dr. Sushil Koirala, BDS, AAAD.



Modification of tooth shape, size & colour

Direct Cosmetic Veneers


Repair of direct and indirect prosthesis

Diastema Closure

## LIFE-LIKE AESTHETICS, UNPARALLELED RESULTS

### Naturally Vivid & Stable Shades

**BEAUTIFIL II LS** is formulated to simulate the optical nature of natural teeth by harmonizing light transmission and diffusion seen in dentin and enamel in any ambient light condition and even reflects natural fluorescence.



Photographs courtesy of Dr. Sushil Koirala, BDS, AAAD.

Photographs courtesy of Dr. Kozo Yoshino, Fukuishi Dental Clinic, Japan.


Photographs courtesy of Dr. Sven Rinke, Hanau, Germany.

The refractive index of the special pre-polymerized and multifunctional fillers is optimized with the resin matrix to match light translucency, hue and brightness of natural teeth.

### Simple & Logical Shade Spectrum


**BEAUTIFIL II LS** offers a practical shade line up that exhibits 3 levels of translucency to easily mimic the enamel as well as deep and superficial dentin layers of natural teeth.

#### BEAUTIFIL II LS






The shades A1-C3 correspond to the VITA® Classical Shades.

### Complementary Line Up BEAUTIFIL II Gingiva



#### BEAUTIFIL II Enamel




### Ideal for Restoration of


- 1) Class I through Class VI cavities
- 2) Modification of tooth shape, size & colour
- 3) Direct Cosmetic Veneers
- 4) Repair of direct and indirect prosthesis
- 5) Diastema Closure

### Exceptional Shade Stability

Less pre and post-cure transition in shade and translucency



Gloss (30 sec.) tested to ISO 2813:2014



### Sustained Surface Luster

Uniform dispersion of filler complex in **BEAUTIFIL II LS** ensures ease of polishing to a high, lasting gloss within a short time.

### EXCEPTIONAL HANDLING & SCULPTABILITY

#### Easy & Simple to adapt

The balance of the multi-filler phase and innovative SRS monomer creates a non-tacky, moderately firm, sculptable paste that is easy to adapt, without any slump. High tolerance to ambient light facilitates longer working time to easily shape the restoration and recreate fine surface details.



#### BEAUTIFIL II LS Syringe 4 gm

Shades available:-

A1, A2, A3, A3.5, A4, B1, B2, B3, C2, C3, A20 (Opacous A2), A30 (Opacous A3), INC (Incisal), BW (Bleach White)

each  
\$59.95

**Buy 5 get 1 FREE**